

EVACUATION PRIMING TIME



$$\text{Time} = 3600 \frac{V}{S} \ln \frac{P_1}{P_2} \quad (\text{seconds})$$

V = total volume M³
 S = capacity of pump M³/Hr
 Ln = the natural logarithm
 P₁ = initial pressure mb abs
 P₂ = final pressure mb abs

PRIMING TIME (SECONDS) PER 100 LITRES OF SYSTEM

Suction Lift M	LV10-5	LV20-5	LV30-5	LV20-20	LV40-20	LVK40-20
2	8	4	3	3	2	3
4	28	9	6	9	5	9
6	104	17	14	16	10	16
8	-	34	30	28	15	28



LIQUID TRANSFER



For information visit: www.liquivac.co.uk or call: 01642 379500

Tomlinson Hall & Co Ltd, Lagonda Road, Cowpen Lane Industrial Estate, Billingham, TS23 4JA

LIQUID TRANSFER



DETAILS OF L282 AVAILABLE ON REQUEST